Instructions for use

Sterile Hypodermic Syringes for Single Use

Please read the instruction book carefully before using 

Date of issue and edition 2025/09 rev1

Product name Sterile hypodermic syringe for single use

Model and specification

Product specification		
Type	Capacity	Luer conical
Three-Piece Type Two-piece type	1ml	Luer lock, Or luer slip
	2ml	
	2.5ml	
	3ml	
	5ml	
	10ml	
	20ml	
	30ml	
	50ml	
	60ml	

Intended use

This device intended to inject fluids into, or withdraw fluids from the body.

Intended patient population

Patients who require injections or withdrawals of medicines and other fluids.

Intended users

Medical professionals.

Indications

None.

Limitation

Use is restricted to trained physicians or caregivers.

Structural composition

The product is composed of outer barrel and plunger (or with piston).

Contra-indication

No contraindications.

Caution

1. Non-toxic, pyrogen-free
2. Destroy if the package is not intact
3. Destroy after single use, harmless treatment shall be carried out in accordance with the regulations on the disposal of medical waste.
4. Not to be used with paraldehyde.

Warning

No warning.

Product tolerance on any graduated capacity and scale interval

Nominal capacity of syringe v ml	Tolerance on any graduated capacity		Maximum scale interval ml
	Less than half nominal capacity	Equal to or greater than half nominal capacity	
$v < 2$	$\pm(1.5\% \text{ of } v + 2\% \text{ of expelled volume})$	$\pm 5\% \text{ of expelled volume}$	0.05
$2 \leq v < 5$	$\pm(1.5\% \text{ of } v + 2\% \text{ of expelled volume})$	$\pm 5\% \text{ of expelled volume}$	0.2
$5 \leq v < 10$	$\pm(1.5\% \text{ of } v + 1\% \text{ of expelled volume})$	$\pm 4\% \text{ of expelled volume}$	0.5
$10 \leq v < 20$	$\pm(1.5\% \text{ of } v + 1\% \text{ of expelled volume})$	$\pm 4\% \text{ of expelled volume}$	1.0
$20 \leq v < 30$	$\pm(1.5\% \text{ of } v + 1\% \text{ of expelled volume})$	$\pm 4\% \text{ of expelled volume}$	2.0
$30 \leq v < 50$	$\pm(1.5\% \text{ of } v + 1\% \text{ of expelled volume})$	$\pm 4\% \text{ of expelled volume}$	2.0
$v \geq 50$	$\pm(1.5\% \text{ of } v + 1\% \text{ of expelled volume})$	$\pm 4\% \text{ of expelled volume}$	5.0

The main performance characteristics device:

1. Dead space: Shall comply with the terms of ISO 7886-1
2. Air tightness (negative pressure):
When tested in accordance with ISO 7886, there shall be no leakage of air past the piston or seal(s), and there shall be no fall in the manometer reading.
3. Water tightness(positive pressure):
When tested in accordance with ISO 7886, there shall be no leakage of water past the piston or seal(s).
4. Fit of plunger stopper/plunger in barrel:
When the syringe is filled with water to the nominal capacity and held vertically with first one end and then the other end uppermost, the piston shall not move by reason of its own mass and the water contained
5. Please refer to ISO 7886-1:2017 for the rest of the performance description.

[Shelf life]: See the label for details.

[Storage and transportation]:

Storage: The sterilized syringe should be stored in a non-corrosive gas, well-ventilated and clean room, and the syringe should be fully protected.

Transportation: The outer packaging should be able to ensure that it will not be damaged under normal transportation conditions.

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Zhejiang Longde Pharmaceutical Co., Ltd.